

# VITROCELL® Smoking Machine VC 1/7

Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		No	Software Option (Bell-shaped Profile) Smoke Box	Conventional Cigarettes	Manually
Electronic cigarettes		Yes	—	Standard E-Cigarettes	Automatically
ENDS cigarette (Electronic Nicotine Delivery)		No	Seals depending on products	Differnt ENDS products possible (e.g. Juul, blu, Alto)	Draw actuated
Renova Zero (ENDS product)		No	—	Renova Zero product	Draw actuated
SMOK (ENDS product)		No	—	SMOK product	Draw actuated
IQOS (Heated Tobacco Product)		No	Software Option (HTP Mode)	IQOS product	Automatically
GLO (Heated Tobacco Product)		No	Software Option (HTP Mode)	GLO product	Automatically
Ploom (Heated Tobacco Product)		No	Software Option (HTP Mode)	Ploom product	Automatically
Vapestarter		No	Software Option (Vapestarter) Rack System	Different HTP products possible	Automatically



VITROCELL® Smoking Machine VC 1/7 standard version



Mode 1: draw actuated cigarettes



Mode 2: vape starters with button activatos for tank products



Mode 3: smoke boxes for combustion cigarettes



Holder system for ENDS products (e.g. blu)