

VITROCELL® Smoking Machine VC 1 S-Type

Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		Yes	—	Conventional Cigarettes	Manually
Electronic cigarettes		No	Software Option (Square Profile)	Standard draw actuated E-Cigarettes	Draw actuated
ENDS cigarette (Electronic Nicotine Delivery)		No	Software Option (Square Profile) Seals depending on product	Differnt ENDS products possible (e.g. Juul, blu, Alto)	Draw actuated
Renova Zero (ENDS product)		No	Software Option (Square Profile)	Renova Zero product	Draw actuated
SMOK (ENDS product)		No	Software Option (Square Profile)	SMOK product	Draw actuated
IQOS (Heated Tobacco Product)		No	Software Option (HTP Mode)	IQOS product	Manually
GLO (Heated Tobacco Product)		No	Software Option (HTP Mode)	GLO product	Manually
Ploom (Heated Tobacco Product)		No	Software Option (HTP Mode)	Ploom product	Manually
Vapestarter		No	Software Option (Vapestarter)	Different HTP products possible	Automatically
Cigar Holder System		No	—	Cigar like products OD 6-22 mm	Manually



VITROCELL® Smoking Machine VC 1 S-Type standard version



Smoke Box and Impinger for Conventional Cigarettes



Option: holder system for ENDS products (e.g. blu)



Option: for ENDS products
(e.g. button or draw actuated and tank products)



Option: rack for ENDS products