


VITROCELL® Smoking Machine VC 1

Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		Yes	Smoke Box	Conventional Cigarettes	Manually
Electronic cigarettes		No	Software Option (Square Profile)	Standard draw actuated E-Cigarettes	Draw actuated
ENDS cigarette (Electronic Nicotine Delivery)		No	Software Option (Square Profile) Seals depending on product	Different ENDS products possible (e.g. Juul, blu, Alto)	Draw actuated
Renova Zero (ENDS product)		No	Software Option (Square Profile)	Renova Zero product	Draw actuated
SMOK (ENDS product)		No	Software Option (Square Profile)	SMOK product	Draw actuated
IQOS (Heated Tobacco Product)		No	Software Option (HTP Mode)	IQOS product	Manually/ Automatically
GLO (Heated Tobacco Product)		No	Software Option (HTP Mode)	GLO product	Manually/ Automatically
Ploom (Heated Tobacco Product)		No	Software Option (HTP Mode)	Ploom product	Manually/ Automatically
Vapestarter		No	Software Option (Vapestarter)	Different HTP products possible	Automatically
Cigar Holder System		No	Smoke Box	Cigar like products OD 6-22 mm	Manually




VITROCELL® Smoking Machine VC 1 Standard version



Standard: for Conventional Cigarettes



Option: for draw actuated e-cigarettes



Option: for Heated Tobacco Products (HTP)



Option: for button actuated e-cigarettes



Option: for any puff-actuated device