


VITROCELL® Smoking Robot VC 10 S-Type


Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		Yes	—	Conventional Cigarettes	Automatically
Electronic cigarettes		No	Software Option (Square Profile) Additional holder	Standard E-Cigarettes	Automatically
ENDS cigarette (Electronic Nicotine Delivery)		No	Software Option (Square Profile) Seals depending on products	Differnt ENDS products possible (e.g. Juul, blu, Alto)	Automatically
Renova Zero (ENDS product)		No	Software Option (Square Profile) Rack System	Renova Zero product	Automatically
SMOK (ENDS product)		No	Software Option (Square Profile) Rack System	SMOK product	Automatically
Heated Tobacco Products (HTP)		No	Including different holder Software Option	Different HTP products	Automatically
IQOS (Heated Tobacco Product)		No	Rack HTP Software	IQOS product	Manually
GLO (Heated Tobacco Product)		No	Rack HTP Software	GLO product	Manually
Ploom (Heated Tobacco Product)		No	Rack System	Ploom product	Manually
Vapestarter		No	Rack System Software Option (Square Profile)	Different HTP products possible	Automatically




VITROCELL® Smoking Robot VC 10 S-Type Standard version




Option: quick-lock holder for HTP products



Option: quick-lock holder for HTP products



Holder for cigarettes



Option: holder for ENDS products



Option: holder system for ENDS products (e.g. blu)