

# VITROCELL® Smoking Robot VC 10

Holder	Product Picture	Standard	Necessary Accessories	Products	Activation
Cigarettes		Yes	—	Conventional Cigarettes	Automatically
Electronic cigarettes		No	Software Option (Square Profile) Additional E-Cig holder	Standard E-Cigarettes	Draw actuated
ENDS cigarette (Electronic Nicotine Delivery)		No	Rack System Seals depending on products	Differnt ENDS products possible (e.g. Juul, blu, Alto)	Draw actuated
Renova Zero (ENDS product)		No	Rack System	Renova Zero product	Draw actuated
SMOK (ENDS product)		No	Rack System	SMOK product	Draw actuated
IQOS (Heated Tobacco Product)		No	Software Option (HTP Mode) Rack System	IQOS product	Manually
GLO (Heated Tobacco Product)		No	Software Option (HTP Mode) Rack System	GLO product	Manually
Ploom (Heated Tobacco Product)		No	Software Option (HTP Mode) Rack System	Ploom product	Manually
Vapestarter		No	Software Option (Vapestarter)	Different HTP products possible	Automatically



VITROCELL® Smoking Robot VC 10 standard version



Option: holder system for ENDS products (e.g. blu)



Standard: automatically lighted cigarettes